

MEDITEST

New technology



MEDITEST 600 / MEDITEST 1300

Froilabo Firlabo


Stable
Reliable
Uniform in T° and Hr

Stability testing climatic cabinet



	MEDITEST 600	MEDITEST 1300
Temperature range	from +15°C up to +60°C	
Controlled relative humidity	from 40 to 90% HR	
actual and setpoint display	2	
Alarm	High and low temperature, high and low humidity	
Security	Independant 3.1 overtemperature safety thermostat	
Capacity (l)	600	1350
External dimensions W x D x H (mm)	700 x 830 x 2130	1430 x 830 x 2130
Internal dimensions W x D x H (mm)	530 x 660 x 1550	1270 x 660 x 1550
Shelves W x D 530 x 655 mm	4	2 x 4
35 mm cable entry port	Right side	
Internal and external cabinet	Stainless steel	
Key locking door	Yes	
Castors	4 (2 with brake)	5 (2 with brake)
CFC and HCFC free	R134a	
Water consumption	0.3 to 0.5 L/24h	
Power supply	230 V - 50 Hz single phase	
Nominal power	1650 W	2400 W

In accordance with ICH guideline Q1A



Long-term condition :	25°C / 60% HR
	30°C / 65% HR
Accelerated condition :	40°C / 75% HR
Intermediate condition :	30°C / 60% HR
Semi permeable packaging	
Long-term condition :	25°C / 40% HR
Accelerated condition :	40°C / 25% HR
Temperature uniformity :	+/-2°C
Humidity uniformity :	+/-5% Hr



8 Bd Monge - 69330 Meyzieu - France Tél: +33 4 78 04 75 75 - Fax: +33 4 78 04 75 76.
export@froilabo.com

www.froilabo.com

This document does not represent a contractual commitment; we reserve the right to make changes without notice