

Electrothermal CMUA Stirring Mantles



The Electrothermal CMUA Stirring Mantles are available in sizes 50 to 5000ml

Electrothermal CMUA Stirring Mantles are designed for stirring and heating of liquid

- Flexible coiled heating element absorb shock, minimizing the risk of flask breaking
- Heating element is attached to a rigid thermal insulating cartridge to ensure maximum heat transfer and support.
- Heating elements and insulation form an easy –to-replace heating cartridge
- Unique air flow through ventilation slots beneath and around the rim and heating cartridge keeps the exterior “cool-to-the-touch”
- Stainless-steel screen covering the heating element is grounded directly to the cable for added protection
- Mantles are double fused for added safety
- Non-skid feet and support clamps add stability

Certifications: CE, CSA

Specifications		
Material		Aluminium case stove painted
Max. Element Temperature		450° C (842° F)
Cat. No.	Capacity	Electrical Requirements
CMUA0050/CE	50ml (0.01gal)	230v 50/60Hz, 60w
CMUA0050/CEX1	50ml (0.01gal)	115v 50/60Hz, 60w
CMUA0050/CEX6	50ml (0.01gal)	230v 50/60Hz, 60w, EU Plug
CMUA0100/CE	100ml (0.02gal)	230v 50/60Hz, 75w
CMUA0100/CEX1	100ml (0.02gal)	115v 50/60Hz, 75w
CMUA0100/CEX6	100ml (0.02gal)	230v 50/60Hz, 75w, EU Plug
CMUA0250/CE	250ml (0.06gal)	230v 50/60Hz, 200w
CMUA0250/CEX1	250ml (0.06gal)	115v 50/60Hz, 200w
CMUA0250/CEX6	250ml (0.06gal)	230v 50/60Hz, 200w, EU Plug
CMUA0500/CE	500ml (0.13gal)	230v 50/60Hz, 280w
CMUA0500/CEX1	500ml (0.13gal)	115v 50/60Hz, 280w
CMUA0500/CEX6	500ml (0.13gal)	230v 50/60Hz, 280w, EU Plug
CMUA1000/CE	1000ml (0.26gal)	230v 50/60Hz, 380w
CMUA1000/CEX1	1000ml (0.26gal)	115v 50/60Hz, 380w
CMUA1000/CEX6	1000ml (0.26gal)	230v 50/60Hz, 380w, EU Plug
CMJA2000/CE	2000ml (0.53gal)	230v 50/60Hz, 500w
CMUA2000/CEX1	2000ml (0.53gal)	115v 50/60Hz, 500w
CMUA2000/CEX6	2000ml (0.53gal)	230v 50/60Hz, 500w, EU Plug
CMUA3000/CE	3000ml (0.79gal)	230v 50/60Hz, 500w
CMUA3000/CEX1	3000ml (0.79gal)	115v 50/60Hz, 500w
CMUA3000/CEX6	3000ml (0.79gal)	230v 50/60Hz, 500w, EU Plug
CMUA5000/CE	5000ml (1.32gal)	230v 50/60Hz, 800w
CMUA5000/CEX1	5000ml (1.32gal)	115v 50/60Hz, 800w
CMUA5000/CEX6	5000ml (1.32gal)	230v 50/60Hz, 800w, EU Plug

