

Electrothermal ABS Top Immersible Stirrers



ABS Top Immersible Stirrers are completely sealed for use in immersion or high humidity applications. Ideal for stirring in water baths, incubators and environmental chambers.

- Control up to six stir pads with one controller
- Stirring capacities from 400ml to 2L
- Case material is moulded ABS
- Speed range from 0 to 150rpm

Includes: Oval stir bar

Warranty: Three years

Certifications: CE

Technical specification

Stirring range	0 to 150 RPM
Temperature range	- 5 to $\pm 50^{\circ}\text{C}$

Ordering information

Model	Load capacity
PS60002	400ml
PS60003	1000ml
PS60004	2000ml