

## General Purpose Low-Profile Stirrer



**Electrothermal General Purpose Low Profile Stirrer is stable and easy to clean.**

Available with single or multi-place stirring positions.

- Uses Stir Trac\* technology-strong magnetic coupling
- Rated for continuous use and sealed to IP65 standards
- Speed Range from 350 to 2000rpm
- Operating temperature 0 to 50 C 0-90% humidity
- Case material is glass moulded resin

**Includes:** 1 stir bar for each stir position

**Warranty:** Three years

**Certifications:** CE

Specifications					
Stirring Range			350-2000 RPM		
Cat. No.	No. of Stirring Positions	Capacity	Overall L x W x H	Weight	Electrical Requirement/ Plug Type
PS60040	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb.)	230V/UK
PS60040X1	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb.)	115V/US
PS60040X6	1	2L	15x15x7cm (6x6x2.75in)	0.8kg (1.8lb)	230V/EU
PS60057	1	24L	27x27x9.6cm (10.5x10.5x3.5 in)	1.9kg (4.2lb)	230V/UK
PS60057X1	1	24L	27x27x9.6cm (10.5x10.5x3.5 in)	1.9kg (4.2lb)	115V/US
PS60057X6	1	24L	27x27x9.6cm (10.5x10.5x3.5 in)	1.9kg (4.2lb)	230V/EU
PS60042	4	1L per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.4kg (5.3lb)	230V/UK
PS60042X1	4	1L per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.4kg (5.3lb)	115V/US
PS60042X6	4	1L per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.4kg (5.3lb)	230V/EU
PS60043	5	400ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.7kg (6lb)	230V/UK
PS60043X1	5	400ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.7kg (6lb)	115V/US
PS60043X6	5	400ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	2.7kg (6lb)	230V/EU
PS60044	9	250ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	3.2kg (7lb)	230V/UK
PS60044X1	9	250ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	3.2kg (7lb)	115V/US
PS60044X6	9	250ml per position	27x27x9.6cm (10.5x10.5x3.5 in)	3.2kg (7lb)	230V/EU