

8 Bd Monge, 69330 Meyzieu, FRANCE

☎ 04 78 04 75 75

Fax 04 78 04 75 76

email : [export@froilabo.com](mailto:export@froilabo.com)



# Emotec - 45°C

Temperature range from -20°C to -45°C  
Electronic controller with digital temperature display

Audible and visible alarm system with recharging  
battery power supply

Dry contact for remote alarm (NO/NF)  
Static cooling system by the shelves

Air cooled condenser with safety pressure switch

Fitted with plastic drawers  
Fitted with 4 castors including 2 with break

Key locking door  
Power supply 230V 50Hz mono – CE compliant



Model	Capacity	External dimensions W*D*H	Drawers	Drawers dimensions W*D*H
EMOTEC 1	217 L	660x918x1 560 mm	5	470x400x170 mm
EMOTEC 2	261 L	660x918x1 760 mm	6	470x400x170 mm
EMOTEC 3	305 L	660x918x1 950 mm	7	470x400x170 mm
EMOTEC 4	350 L	660x918x2 150 mm	8	470x400x170 mm
EMOTEC 5	520 L	753x960x2 020 mm	14	290x500x170 mm