

## 13 $\mu$ Mini Fuge

- > Compact design to minimize bench space used
- > Rapid acceleration (<12 sec) and deceleration (<16 sec)
- > Simple display with mono LCD
- > Safe operation with lid-lock function
- > Autoclavable
- > Corrosion resistant rotor
- > Two rotors included: 12 x 1.5/2ml Tubes & 32x0.2ml Tubes or Strips
- > CE Marked



### Technical Specifications

Maximum RPM	13,500 RPM
Maximum RCF	12,300 xg
Maximum Capacity	1.5ml / 2.0ml x 12
Run Time	≤ 30 min
Noise Level	≤ 60dB
Acceleration Time	≤ 12 seconds
Deceleration Time	≤ 16 seconds
Safety	Lid-lock
Available Rotors	1.5ml x 12, 0.2ml x 32 PCR Strips
Power & Frequency	220V, 50/60 Hz
Dimensions (WxDxH) mm	208 x 245 x 145
Weight	4.4kg
CE Marking	Yes