

Electrothermal Paraffin Wax Dispenser

The Electrothermal Paraffin Wax Dispenser provides immediate “Wax On Tap” for Pathology laboratories etc.



The problem of wax solidifying at the tap is completely eliminated by the incorporation of a thermally insulated tap heater with independent temperature control.

Capacity:

- 4.5 litres (1 gallon)

Temperature:

- Ambient to 75°C.
- Fitted with safety over-temperature cut-out at 105°C and non-resettable thermal fuse.

Temperature Control:

- A built-in thermostat provides temperature control from ambient 75°C.
- Tap Heater maintains temperature within the tap, keeping the wax molten.

Tap:

- Swivel lever non-drip with replaceable filter

Warranty: One year, parts and labour

Certifications: CE

Specifications		
Material	Stove Painted Aluminium Case	
Max. Element Temperature	5°C to 40°C. Humidity not to exceed 80%.	
Cat. No.	Capacity	Electrical Requirement
MH8523B	4.5 litres (1 gallon)	230V 50/60Hz
MH8523BX1	4.5 litres (1 gallon)	115v 50/60Hz
MH8523BX6	4.5 litres (1 gallon)	230v 50/60Hz, EU Plug